CBDepot, s.r.o., Masarykova třída 1595/54 CZ-415 01 Teplice, Czech Republic	Doc. Number: PSCH1	CBDepot.eu Upmarket Carmabinoid Solutions
Purchasing Specification Cannabis Herba	Issue Number: 13	
Prepared by: Michal Tőzsér	Created on: 31/03/22	

PURCHASING SPECIFICATION

Cannabis Herba

the principles 2 95x95cm and tacked on a 1 ed lab on delivation and tacked on a 1 ed lab on delivation and tacked and securited lab on delivation and tacked lab	handling and processing, with traceability of the estated in EMA-GACP (Good Agricultural and destated in EMA-G
the principles 2 95x95cm and tacked on a 1 ed lab on delivation and tacked on a 1 ed lab on delivation and tacked and securited lab on delivation and tacked lab	se stated in EMA-GACP (Good Agricultural and d height 150cm without palette. Big bag fit with 4 120x80cm or 120x100cm palette Each big bag marked ivery required. RE, each bag marked with net weight, tied or stitched, red. Maximum height of the palette 170cm. Each livery required. Method Visual Ph.Eur. 2.8.2 In-House GC FID Total ash (Ph.Eur.2.4.16) Los on drying (Ph.Eur.2.2.32) Drying at max. 60°C In-House GC FID
the principles 2 95x95cm and tacked on a 1 ed lab on delivation and tacked on a 1 ed lab on delivation and tacked and securited lab on delivation and tacked lab	se stated in EMA-GACP (Good Agricultural and d height 150cm without palette. Big bag fit with 4 120x80cm or 120x100cm palette Each big bag marked ivery required. RE, each bag marked with net weight, tied or stitched, red. Maximum height of the palette 170cm. Each livery required. Method Visual Ph.Eur. 2.8.2 In-House GC FID Total ash (Ph.Eur.2.4.16) Los on drying (Ph.Eur.2.2.32) Drying at max. 60°C In-House GC FID
tacked on a 1 ed lab on delivith same TAF ped and secur ted lab on delivith same TAF ped and secur ted lab on delivith same TAF ped and secur ted lab on delivith same TAF ped and secur ted lab on delivith same secure ted l	120x80cm or 120x100cm palette Each big bag marked ivery required. RE, each bag marked with net weight, tied or stitched, red. Maximum height of the palette 170cm. Each livery required. Method Visual Ph.Eur. 2.8.2 In-House GC FID Total ash (Ph.Eur.2.4.16) Los on drying (Ph.Eur.2.2.32) Drying at max. 60°C In-House GC FID In-House GC FID In-House GC FID In-House GC FID
tacked on a 1 ed lab on delivith same TAF ped and secur ted lab on delivith same TAF ped and secur ted lab on delivith same TAF ped and secur ted lab on delivith same TAF ped and secur ted lab on delivith same secure ted l	120x80cm or 120x100cm palette Each big bag marked ivery required. RE, each bag marked with net weight, tied or stitched, red. Maximum height of the palette 170cm. Each livery required. Method Visual Ph.Eur. 2.8.2 In-House GC FID Total ash (Ph.Eur.2.4.16) Los on drying (Ph.Eur.2.2.32) Drying at max. 60°C In-House GC FID In-House GC FID In-House GC FID In-House GC FID
- [% w/w] - [% w/w]	Visual Ph.Eur. 2.8.2 In-House GC FID Total ash (Ph.Eur.2.4.16) Los on drying (Ph.Eur.2.2.32) Drying at max. 60°C In-House GC FID
- [% w/w] - [% w/w]	Visual Ph.Eur. 2.8.2 In-House GC FID Total ash (Ph.Eur.2.4.16) Los on drying (Ph.Eur.2.2.32) Drying at max. 60°C In-House GC FID
[% w/w] - [% w/w]	Ph.Eur. 2.8.2 In-House GC FID Total ash (Ph.Eur.2.4.16) Los on drying (Ph.Eur.2.2.32) Drying at max. 60°C In-House GC FID
[% w/w]	In-House GC FID Total ash (Ph.Eur.2.4.16) Los on drying (Ph.Eur.2.2.32) Drying at max. 60°C In-House GC FID
[% w/w]	Total ash (Ph.Eur.2.4.16) Los on drying (Ph.Eur.2.2.32) Drying at max. 60°C In-House GC FID
[% w/w]	Los on drying (Ph.Eur.2.2.32) Drying at max. 60°C In-House GC FID
[% w/w] [% w/w] [% w/w] [% w/w] [% w/w] [% w/w]	In-House GC FID
[% w/w] [% w/w] [% w/w] [% w/w] [% w/w]	In-House GC FID In-House GC FID In-House GC FID In-House GC FID
[% w/w] [% w/w] [% w/w] [% w/w] [% w/w]	In-House GC FID In-House GC FID In-House GC FID In-House GC FID
[% w/w] [% w/w] [% w/w] [% w/w]	In-House GC FID In-House GC FID In-House GC FID
[% w/w] [% w/w] [% w/w]	In-House GC FID In-House GC FID
[% w/w] [% w/w]	In-House GC FID
[% w/w] [% w/w]	In-House GC FID
[% w/w]	
ppm	In-House GC FID
	Ph. Eur. 2.4.27
ppm	Ph. Eur. 2.4.27
ppm	Ph. Eur. 2.4.27
ppm	Ph. Eur. 2.4.27
CFU/g	
CFU/g	Microbiological quality of herbal medicinal product for oral use and extracts used in their preparation
CFU/g	(Ph. Eur.5.1.8,
-	2.6.12, 2.6.13, 2.6.31)
ppb	Ph. Eur. 2.8.18
ppb	Ph. Eur 2.8.22
ppm	Pesticide residues (Ph. Eur.2.8.13)
	(i ii. Edi. 2.0.13)
ned for Seller me, position nature, stamp	
r	ppb ppm ned for Selle